
PyAbacus Documentation

Release 1.1.0

Tausand Electronica

Apr 27, 2020

CONTENTS

1	Contents	3
1.1	pyAbacus.core	3
1.2	pyAbacus.exceptions	7
1.3	pyAbacus.constants	8
2	Indices and tables	11
	Python Module Index	13
	Index	15



pyAbacus was built to simplify the usage of **Tausand** Abacus family of coincidence counters, providing a library aimed to interface these devices using Python coding.

CONTENTS

1.1 pyAbacus.core

```
class pyAbacus.core.AbacusSerial (port)
    Builds a serial port from pyserial.

    findIdn ()
        Requests the device for its string identifier (IDN) using serial port.

    flush ()

    getIdn ()
        Gets the device string identifier (IDN) from local memory.

    getNChannels ()
        Gets the number of input channels in the device.

    readSerial ()

    testDevice ()

    writeSerial (command, address, data_16o32)

class pyAbacus.core.CountersValues (n_channels)
    Keeps a set of measurements from counters within a device.

    getCountersID ()
        Gets the counters_id (consecutive number of measurements) field from a set of measurements.

    getNumericAddresses ()

    getTimeLeft ()
        Gets the time_left (time in ms for next measurement to be available) field from a set of measurements.

    getValue (channel)

    getValues (channels)

    getValuesFormatted (channels)

    setCountersID (id)

    setTimeLeft (time)

    setValueFromArray (address, value)

    time_left = None
        in ms

class pyAbacus.core.Settings2Ch
```

getAddressAndValue (*timer*)

getSetting (*timer*)

getSettingStr (*timer*)

setSetting (*setting, value*)

class `pyAbacus.core.Settings48Ch`

4 and 8 channel devices use as time base a second. Nevertheless 2 channel uses ns for all timers with the exception of the sampling time (ms).

exponentRepresentationToValue (*c, e*)

exponentsToBits (*c, e*)

fromBitsToValue (*bits*)

getAddressAndValue (*timer*)

getChannels ()

getSetting (*timer*)

For all timers: returns nanoseconds, for sampling returns ms.

getSettingStr (*timer*)

initAddresses ()

setSetting (*setting, value*)

For all timers: value is in nanoseconds, for sampling in ms.

valueToExponentRepresentation (*number*)

class `pyAbacus.core.Settings4Ch`

4 and 8 channel devices use as time base a second. Nevertheless 2 channel uses ns for all timers with the exception of the sampling time (ms).

class `pyAbacus.core.Settings8Ch`

4 and 8 channel devices use as time base a second. Nevertheless 2 channel uses ns for all timers with the exception of the sampling time (ms).

class `pyAbacus.core.Stream` (*abacus_port, counters, output_function=<built-in function print>*)

setCounters (*counters*)

start ()

stop ()

`pyAbacus.core.close` (*abacus_port*)

Closes a Tausand Abacus device session

`pyAbacus.core.dataArraysToCounters` (*abacus_port, addresses, data*)

Saves in local memory the values of device's counters.

Args: *abacus_port*: device port.

addresses: list of integers with device's register addresses.

data: list of integers with device's register values.

Returns: List of counter values as registered within the device.

`pyAbacus.core.dataArraysToSettings` (*abacus_port, addresses, data*)

Saves in local memory the values of device's settings.

Args: abacus_port: device port.

addresses: list of integers with device's register addresses.

data: list of integers with device's register values.

Returns: List of settings as registered within the device.

`pyAbacus.core.dataStreamToDataArrays (input_string, chunk_size=3)`

Builds data from string read on serial port.

Args: input_string: stream of bytes to convert. Should have the appropriate format, as given by a Tausand Abacus device.

chunk_size : integer, number of bytes per single data row.

- Use chunk_size=3 for devices with inner 16-bit registers e.g. Tausand Abacus AB1002, where byte streams are: {address,MSB,LSB}.

- Use chunk_size=5 for devices with inner 32-bit registers e.g. Tausand Abacus AB1004, where byte streams are: {address,MSB,2nd-MSB,2nd-LSB,LSB}.

Returns: Two lists of integer values: addresses, data.

Raises: AbacusError: Input string is not valid chunk size must either be 3 or 5.

`pyAbacus.core.findDevices (print_on=True)`

Returns a list of connected and available devices that match with a Tausand Abacus.

Scans all serial ports, and asks each of them their descriptions. When a device responds with a valid string, e.g. "Tausand Abacus AB1002", the port is included in the final answer.

Args: print_on: bool When True, prints devices information.

Returns: ports, len(ports) List of valid ports, and its length.

`pyAbacus.core.getAllCounters (abacus_port)`

Reads all counters from a Tausand Abacus device.

With a single call, this function reads all the counters within the device, including single-channel counters, 2-fold coincidence counters and multi-fold coincidence counters.

Example:

```
counters, counters_id = getAllCounters('COM3')
```

Reads data from the device in port 'COM3', and might return for example,

```
counters = {A:1023, B:1038, AB: 201}
```

```
counters.A = 1023
```

```
counters.B = 1038
```

```
counters.AB = 201
```

```
counters_id = 37
```

meaning that this is the 37th measurement made by the device, and the measurements were 1023 counts in A, 1038 counts in B, and 201 coincidences between A and B.

Args: abacus_port: device port.

Returns: CountersValues class object including counter values as registered within the device, and the sequential number of the reading.

`pyAbacus.core.getAllSettings(abacus_port)`

Reads all settings from a Tausand Abacus device.

With a single call, this function reads all the settings within the device, including sampling time, coincidence window, delay per channel and sleep time per channel.

Example: `settings = getAllSettings('COM3')`

Reads settings from the device in port 'COM3', and might return for example,

```
delay_A (ns): 0
delay_B (ns): 20
sleep_A (ns): 0
sleep_B (ns): 0
coincidence_window (ns): 10
sampling (ms): 1300
```

Args: `abacus_port`: device port.

Returns: List of settings as registered within the device.

`pyAbacus.core.getChannelsFromName(name)`

Returns the number of input channels by reading the device name.

For example, if `name="Tausand Abacus AB1004"`, returns 4.

Args: `name`: idn string of the device.

Returns: integer, number of input channels in device.

Raises: `AbacusError`: Not a valid abacus.

`pyAbacus.core.getCountersID(abacus_port)`

Reads the *counters_id* (consecutive number of measurements) in a Tausand Abacus.

When a new configuration is set, *counters_id*=0, indicating no valid data is available.

Each time a new set of valid measurements is available, *counters_id* increments 1 unit.

counters_id overflows at 1 million, starting over at *counters_id*=1.

Args: `abacus_port`: device port.

Returns: integer, *counters_id* value.

`pyAbacus.core.getFollowingCounters(abacus_port, counters)`

`pyAbacus.core.getIdn(abacus_port)`

Reads the identifier string model (IDN) from a Tausand Abacus.

Example: `myidn = getIdn('COM3')`

might return

```
myidn = "Tausand Abacus AB1002"
```

Args: abacus_port: device port.

Returns: IDN string.

`pyAbacus.core.getSetting(abacus_port, setting)`

Get a single configuration setting within a Tausand Abacus.

Args: abacus_port: device port

setting: name of the setting to be written. Valid strings are: “sampling”, “coincidence_window”, “delay_N”, “sleep_N”, where “N” refers to a channel (A,B,C,D,...).

Returns: value for the setting. For “sampling”, value in ms; for other settings, value in ns.

`pyAbacus.core.getTimeLeft(abacus_port)`

Reads the remaining time for the next measurement to be ready, in ms.

Args: abacus_port: device port

Returns: integer, in ms, of time left for next measurement.

`pyAbacus.core.open(abacus_port)`

Opens a session to a Tausand Abacus device

`pyAbacus.core.readSerial(abacus_port)`

Reads bytes available at the specified serial port.

`pyAbacus.core.renameDuplicates(old)`

`pyAbacus.core.setAllSettings(abacus_port, new_settings)`

`pyAbacus.core.setSetting(abacus_port, setting, value)`

Sets a configuration setting within a Tausand Abacus.

Example: `setSetting('COM3', 'sampling', 1300)`

sets the sampling time to 1300 ms to a device in port ‘COM3’.

Args: abacus_port: device port

setting: name of the setting to be written. Valid strings are: “sampling”, “coincidence_window”, “delay_N”, “sleep_N”, where “N” refers to a channel (A,B,C,D,...).

value: new value for the setting. For “sampling”, value in ms; for other settings, value in ns.

`pyAbacus.core.writeSerial(abacus_port, command, address, data_16o32)`

Low level function. Writes in the specified serial port an instruction built based on command, memory address and data.

1.2 pyAbacus.exceptions

exception `pyAbacus.exceptions.AbacusError(message=)`

An unexpected error occurred.

exception `pyAbacus.exceptions.BaseError(message)`

exception `pyAbacus.exceptions.CheckSumError`

An error occurred while doing check sum.

exception `pyAbacus.exceptions.InvalidValueError(message=)`

The selected value is not valid

exception `pyAbacus.exceptions.TimeoutError(message=)`

A time out error occurred

1.3 pyAbacus.constants

`pyAbacus.constants.ADDRESS_DIRECTORY_2CH = {'coincidence_window_ms': 22, 'coincidence_win`
Memory addresses

`pyAbacus.constants.BAUDRATE = 115200`
Default baudrate for the serial port communication

`pyAbacus.constants.BOUNCE_TIMEOUT = 1`
Number of times a specific transmission is tried

`pyAbacus.constants.COINCIDENCE_WINDOW_DEFAULT_VALUE = 10`
Default coincidence window time value (ns).

`pyAbacus.constants.COINCIDENCE_WINDOW_MAXIMUM_VALUE = 10000`
Maximum coincidence window time value (ns).

`pyAbacus.constants.COINCIDENCE_WINDOW_MINIMUM_VALUE = 5`
Minimum coincidence window time value (ns).

`pyAbacus.constants.COINCIDENCE_WINDOW_STEP_VALUE = 5`
Increase ratio on the coincidence window time value (ns).

`pyAbacus.constants.COUNTERS_VALUES = {}`
Global counters values variable

`pyAbacus.constants.CURRENT_OS = 'linux'`
Current operative system

`pyAbacus.constants.DELAY_DEFAULT_VALUE = 0`
Default delay time value (ns).

`pyAbacus.constants.DELAY_MAXIMUM_VALUE = 100`
Maximum delay time value (ns).

`pyAbacus.constants.DELAY_MINIMUM_VALUE = 0`
Minimum delay time value (ns).

`pyAbacus.constants.DELAY_STEP_VALUE = 5`
Increase ratio on the delay time value (ns).

`pyAbacus.constants.END_COMMUNICATION = 4`
End of message

`pyAbacus.constants.MAXIMUM_WRITING_TRIES = 20`
Number of tries done to write a value

`pyAbacus.constants.READ_VALUE = 14`
Reading operation signal

`pyAbacus.constants.SAMPLING_DEFAULT_VALUE = 1000`
Default sampling time value (ms)

`pyAbacus.constants.SAMPLING_VALUES = [1, 2, 5, 10, 20, 50, 100, 200, 500, 1000, 2000, 5000,`
From (1, 2, 5) ms to 1000 s

`pyAbacus.constants.SETTINGS = {}`
Global settings variable

`pyAbacus.constants.SLEEP_DEFAULT_VALUE = 0`
Default sleep time value (ns).

`pyAbacus.constants.SLEEP_MAXIMUM_VALUE = 100`
Maximum sleep time value (ns).

`pyAbacus.constants.SLEEP_MINIMUM_VALUE = 0`
Minimum sleep time value (ns).

`pyAbacus.constants.SLEEP_STEP_VALUE = 5`
Increase ratio on the sleep time value (ns).

`pyAbacus.constants.START_COMMUNICATION = 2`
Begin message signal

`pyAbacus.constants.TIMEOUT = 0.5`
Maximum time without answer from the serial port

`pyAbacus.constants.WRITE_VALUE = 15`
Writing operation signal

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

p

`pyAbacus.constants`, 8
`pyAbacus.core`, 3
`pyAbacus.exceptions`, 7

A

AbacusError, 7
 AbacusSerial (class in *pyAbacus.core*), 3
 ADDRESS_DIRECTORY_2CH (in module *pyAbacus.constants*), 8

B

BaseError, 7
 BAUDRATE (in module *pyAbacus.constants*), 8
 BOUNCE_TIMEOUT (in module *pyAbacus.constants*), 8

C

ChecksumError, 7
 close() (in module *pyAbacus.core*), 4
 COINCIDENCE_WINDOW_DEFAULT_VALUE (in module *pyAbacus.constants*), 8
 COINCIDENCE_WINDOW_MAXIMUM_VALUE (in module *pyAbacus.constants*), 8
 COINCIDENCE_WINDOW_MINIMUM_VALUE (in module *pyAbacus.constants*), 8
 COINCIDENCE_WINDOW_STEP_VALUE (in module *pyAbacus.constants*), 8
 COUNTERS_VALUES (in module *pyAbacus.constants*), 8
 CountersValues (class in *pyAbacus.core*), 3
 CURRENT_OS (in module *pyAbacus.constants*), 8

D

dataArraysToCounters() (in module *pyAbacus.core*), 4
 dataArraysToSettings() (in module *pyAbacus.core*), 4
 dataStreamToDataArrays() (in module *pyAbacus.core*), 5
 DELAY_DEFAULT_VALUE (in module *pyAbacus.constants*), 8
 DELAY_MAXIMUM_VALUE (in module *pyAbacus.constants*), 8
 DELAY_MINIMUM_VALUE (in module *pyAbacus.constants*), 8
 DELAY_STEP_VALUE (in module *pyAbacus.constants*), 8

E

END_COMMUNICATION (in module *pyAbacus.constants*), 8
 exponentRepresentationToValue() (*pyAbacus.core.Settings48Ch* method), 4
 exponentsToBits() (*pyAbacus.core.Settings48Ch* method), 4

F

findDevices() (in module *pyAbacus.core*), 5
 findIdn() (*pyAbacus.core.AbacusSerial* method), 3
 flush() (*pyAbacus.core.AbacusSerial* method), 3
 fromBitsToValue() (*pyAbacus.core.Settings48Ch* method), 4

G

getAddressAndValue() (*pyAbacus.core.Settings2Ch* method), 3
 getAddressAndValue() (*pyAbacus.core.Settings48Ch* method), 4
 getAllCounters() (in module *pyAbacus.core*), 5
 getAllSettings() (in module *pyAbacus.core*), 6
 getChannels() (*pyAbacus.core.Settings48Ch* method), 4
 getChannelsFromName() (in module *pyAbacus.core*), 6
 getCountersID() (in module *pyAbacus.core*), 6
 getCountersID() (*pyAbacus.core.CountersValues* method), 3
 getFollowingCounters() (in module *pyAbacus.core*), 6
 getIdn() (in module *pyAbacus.core*), 6
 getIdn() (*pyAbacus.core.AbacusSerial* method), 3
 getNChannels() (*pyAbacus.core.AbacusSerial* method), 3
 getNumericAddresses() (*pyAbacus.core.CountersValues* method), 3
 getSetting() (in module *pyAbacus.core*), 7
 getSetting() (*pyAbacus.core.Settings2Ch* method), 4
 getSetting() (*pyAbacus.core.Settings48Ch* method), 4

`getSettingStr()` (*pyAbacus.core.Settings2Ch method*), 4
`getSettingStr()` (*pyAbacus.core.Settings48Ch method*), 4
`getTimeLeft()` (*in module pyAbacus.core*), 7
`getTimeLeft()` (*pyAbacus.core.CountersValues method*), 3
`getValue()` (*pyAbacus.core.CountersValues method*), 3
`getValues()` (*pyAbacus.core.CountersValues method*), 3
`getValuesFormatted()` (*pyAbacus.core.CountersValues method*), 3

I

`initAddreses()` (*pyAbacus.core.Settings48Ch method*), 4
`InvalidValueError`, 7

M

`MAXIMUM_WRITING_TRIES` (*in module pyAbacus.constants*), 8
`module`
 pyAbacus.constants, 8
 pyAbacus.core, 3
 pyAbacus.exceptions, 7

O

`open()` (*in module pyAbacus.core*), 7

P

`pyAbacus.constants`
 module, 8
`pyAbacus.core`
 module, 3
`pyAbacus.exceptions`
 module, 7

R

`READ_VALUE` (*in module pyAbacus.constants*), 8
`readSerial()` (*in module pyAbacus.core*), 7
`readSerial()` (*pyAbacus.core.AbacusSerial method*), 3
`renameDuplicates()` (*in module pyAbacus.core*), 7

S

`SAMPLING_DEFAULT_VALUE` (*in module pyAbacus.constants*), 8
`SAMPLING_VALUES` (*in module pyAbacus.constants*), 8
`setAllSettings()` (*in module pyAbacus.core*), 7
`setCounters()` (*pyAbacus.core.Stream method*), 4
`setCountersID()` (*pyAbacus.core.CountersValues method*), 3

`setSetting()` (*in module pyAbacus.core*), 7
`setSetting()` (*pyAbacus.core.Settings2Ch method*), 4
`setSetting()` (*pyAbacus.core.Settings48Ch method*), 4
`setTimeLeft()` (*pyAbacus.core.CountersValues method*), 3
`SETTINGS` (*in module pyAbacus.constants*), 8
`Settings2Ch` (*class in pyAbacus.core*), 3
`Settings48Ch` (*class in pyAbacus.core*), 4
`Settings4Ch` (*class in pyAbacus.core*), 4
`Settings8Ch` (*class in pyAbacus.core*), 4
`setValueFromArray()` (*pyAbacus.core.CountersValues method*), 3
`SLEEP_DEFAULT_VALUE` (*in module pyAbacus.constants*), 8
`SLEEP_MAXIMUM_VALUE` (*in module pyAbacus.constants*), 8
`SLEEP_MINIMUM_VALUE` (*in module pyAbacus.constants*), 9
`SLEEP_STEP_VALUE` (*in module pyAbacus.constants*), 9
`start()` (*pyAbacus.core.Stream method*), 4
`START_COMMUNICATION` (*in module pyAbacus.constants*), 9
`stop()` (*pyAbacus.core.Stream method*), 4
`Stream` (*class in pyAbacus.core*), 4

T

`testDevice()` (*pyAbacus.core.AbacusSerial method*), 3
`time_left` (*pyAbacus.core.CountersValues attribute*), 3
`TIMEOUT` (*in module pyAbacus.constants*), 9
`TimeOutError`, 7

V

`valueToExponentRepresentation()` (*pyAbacus.core.Settings48Ch method*), 4

W

`WRITE_VALUE` (*in module pyAbacus.constants*), 9
`writeSerial()` (*in module pyAbacus.core*), 7
`writeSerial()` (*pyAbacus.core.AbacusSerial method*), 3